

Packing List

Reference: SC-018-006

Shipper / Exporter:

Sercon Limited
Unit 3B Crewe Trade park
Gateway
Crewe
CW1 6JT
UK

Delivery address :

E.S. Tech KYOTO
Research Institute for Production Development 8th
Institute
15 Morimoto-cho, Shimogamo, Sakyo-ku
Kyoto, 606-0805 Japan
TEL : 81 (80) 3036-5088

Consignee :

Taiyo Nippon Sanso Corporation
Stable Isotope Div.
Fujj Bldg. 5-30-9 Shiba, Minato-Ku
Tokyo 108-0014 Japan
Tel : 81(3)-5439-5897

Tel: +44 (0) 1270 580008
Reference: IN1268

Contact Name: Yasuo Ogawa
Contact Tel Number: +81 (3) 5439 5897

BOX: 1

Quantity	Product Code	Product Description	Serial Number: if available
☒ 1	20602/1	ABCA	20602/2-2021/1268

Packed by:

James Haver

Date:

9.12.21

Claims for damage or shortages must be made within 3 days

Packing List

Reference: SC-018-006

Shipper / Exporter:

Sercon Limited
 Unit 3B Crewe Trade park
 Gateway
 Crewe
 CW1 6JT
 UK

Delivery address :

E.S. Tech KYOTO
 Research Institute for Production Development 8th
 Institute
 15 Morimoto-cho, Shimogamo, Sakyo-ku
 Kyoto, 606-0805 Japan
 TEL : 81 (80) 3036-5088

Consignee :

Taiyo Nippon Sanso Corporation
 Stable Isotope Div.
 Fuji Bldg. 5-30-9 Shiba, Minato-Ku
 Tokyo 108-0014 Japan
 Tel : 81(3)-5439-5897

Tel: +44 (0) 1270 580008
 Reference: IN1268

Contact Name: Yasuo Ogawa
 Contact Tel Number: +81 (3) 5439 5897

BOX: 2

Quantity	Product Code	Product Description	Serial Number: if available
☒ 5		CETAC racks	
☒ 1		CETAC base	
☒ 1		Rotary Pump	200259990
☒ 1		Permanent Magnet + feet	NO 4
☒ 1		Japan Mains Cable (For ABCA) IN CABLE BAG	
☒ 1		Ribbon Cable - to gas box IN CABLE BAG	
☒ 1		Flow Cable - to gas box IN CABLE BAG	
☒ 1	70101-1	ABCA Spares Kit	
☒ 1		Manuals On Software Drive. (In PC box)	
☒ 1		Software Flash Drive (In PC Box)	
☒ 1		Windows Recovery Drive (In PC Box)	
☒ 1		Monitor (1 x Japan mains cable in box)	JBD5WB3
☒ 1		PC (1 x Japan mains cable in box)	CT7OMF3
☒ 3		Connection Tubing	

Packed by:



Date:

9.12.21

Claims for damage or shortages must be made within 3 days

Packing List

Reference: SC-018-006

Shipper / Exporter:

Sercon Limited
Unit 3B Crewe Trade park
Gateway
Crewe
CW1 6JT
UK

Delivery address :

E.S. Tech KYOTO
Research Institute for Production Development 8th
Institute
15 Morimoto-cho, Shimogamo, Sakyo-ku
Kyoto, 606-0805 Japan
TEL : 81 (80) 3036-5088

Consignee :

Taiyo Nippon Sanso Corporation
Stable Isotope Div.
Fuji Bldg. 5-30-9 Shiba, Minato-Ku
Tokyo 108-0014 Japan
Tel : 81(3)-5439-5897

Tel: +44 (0) 1270 580008
Reference: IN1268

Contact Name: Yasuo Ogawa
Contact Tel Number: +81 (3) 5439 5897

BOX: 3

Quantity	Product Code	Product Description	Serial Number: if available
☒ 1		CETAC Autosampler	052114A7400HR
☒ 1		CETAC Arm	
☒ 1		CETAC Needle (In top of box)	
☒ 1		Gas Box	10602/2021/1268
☒ 1		Power supply - CETAC cable	
☒ 3		Mains leads - CETAC cable (Including Japan lead)	
☒ 1		USB Cable - CETAC cable	

Packed by:

James Hansen

Date:

9.12.21

Claims for damage or shortages must be made within 3 days